



www.mesamed.de

ECG cable for Physio Control Lifepak 12, Lifepak 15, Lifepak 20E Defibrillator



P/N MON-120110

Yoke cable type Physio Control Lifepak 12, Lifepak 15, Lifepak 20E Defibrillator, IEC, yoke for 4-lead connection and yoke for 6-lead connection, with resistor, length 2,4 m, Ref.No. 11110-000110 (or Ref.No.11110-000111 with length 1,5 m)

P/N ADD-120110.4

ECG adapter cable type Physio Control Lifepak 12, Lifepak 15, Lifepak 20E Defibrillator,
IEC, with yoke and 4-lead limb cable with snaps, Ref.No. 11110-000104,
for connection to trunk cable P/N MON-120110, Ref.No. 11110-000110

P/N ADD-120110.6

ECG adapter cable type Physio Control Lifepak 12, Lifepak 15, Lifepak 20E Defibrillator,
IEC, with yoke and 6-lead chest cable with snaps, Ref.No. 11110-000105,
for connection to trunk cable P/N MON-120110, Ref.No. 11110-000110

Minimum order qty. is 1 set of MON-120110, ADD-120110.4 and ADD-120110.6


Information:

The trunkcable P/N MON-120110 (Ref.Nr. 11110-000110)
and the connected yoke with 4-lead limb cable with snaps P/N ADD-120110.4
(Ref.Nr. 11110-000104) are the same as the original Lifepak
trunkcable with integrated 4-lead limb cable Ref.Nr. 11111-000021




www.mesamed.de

Compatible ECG Cable – Physio Control Lifepak 12, Lifepak 15, Lifepak 20E Defibrillator



Art.Nr. MON-120110.151

Complete ECG cable type Physio Control Lifepak 12
Lifepak 15, Lifepak 20E Defibrillator, IEC,
with integrated 3-lead limb cable with snap
Ref.No.11110-000030



Art.Nr. MON-120110.251

Complete ECG cable type Physio Control Lifepak 12
Lifepak 15, Lifepak 20E Defibrillator, IEC,
with integrated 4-lead limb cable and
1-lead chest cable, with snap
Ref.No.11110-000067

Lifepak is a registered trademark of Medtronic

MESA Medizintechnik GmbH

Schaerflmuehlweg 4

D-83671 Benediktbeuern

GERMANY

Tel.: +49 (0) 8857 – 6918-0

Fax: +49 (0) 8857 – 6918-29

Email: info@mesamed.de

www.mesamed.de