

Flow Sensors – More Choice

The Flow Sensors range offers you from “**One Hand**” products of a **proven performance**, **guaranteed compatibility at affordable prices**.

- Proven quality
- Versatility in application
- Affordable prices

Proven quality

The rigorous selection of **high quality** medical grade materials used for the production of the Flow Sensors and EN ISO based internal quality system makes it possible to meet and exceed the demands of medical practice for durability.

The unique high quality of our Flow Sensors make this product range **highly appreciated** by our customers world wide.

Versatility in application

The SpiroQuant A+ and SpiroQuant H were designed for a wide range of medical applications in intensive and emergency care as well as in ambulant care.

The EnviteC Flow Sensors range delivers top product that combine high precision, high linearity with small airway resistance. On request we can provide you your own OEM based Flow Sensor.



Affordable prices

Now you can **leave your logistic problems behind you** . We provide the Flow Sensors you need, uncomplicated from one hand.

The affordable price performance will support your budget.



Flow Sensors Cross Reference

Description	Manufactures		Measuring Principle	EnviteC		Application
	Monitors	Sensor PN		Part Number	Connector	
 <p>SpiroQuant A+</p>	Dräger Evita Range	8403735	Constant –temperature-hot wire-anaemometry	1000470	Dräger specific	Compatible with Dräger Spirolog-sensor PN 8403735
 <p>SpiroQuant H</p>	HamiltonH72 VR-1, Galileo, Amadeus, Veolar, Raphael.	H7400 AccuTach™	Differential pressure	07-00-0001	Sensor Airway OD 15mm med cone according (EN-ISO 5356-1) Airway Adapter OD 22mm ID 15 mm med cone according (EN-ISO 5356-1) Tube ID 3mm	Compatible with H7400 AccuTach™ applicable in Hamilton Ventilators; H72 VR-1, Galileo, Amadeus, Veolar, Raphael.