




























Temperature Sensors Serie 400 compatible

Monitor-Connection	Monitor	Part-Number	Patient	Application	Reusable	Disposable	Cable-Length [cm]	Diameter-Cable [mm]	Connector	Compatibility	Diameter-Head [mm]	Sensor Head
	Hewlett Packard, Phillips	HP-AR-275A	Adult	Rectal/Esophageal	x		275	4,0	Hp	YSI 400 Serie	5	
		HP-PR-275A	Pediatric	Rectal/Esophageal	x		275	3,2	Hp	YSI 400 Serie	4	
		HP-NR-275A	Infant/Neonatal	Rectal/Esophageal	x		275	2,3	Hp	YSI 400 Serie	3	
	Hewlett Packard, Phillips	HP-SAR-275A	Adult	Skin-Contact	x		275	4,0	Hp	YSI 400 Serie	18	
		HP-SPR-275A	Pediatric	Skin-Contact	x		275	3,2	Hp	YSI 400 Serie	10	
		HP-SNR-275A	Infant/Neonatal	Skin-Contact	x		275	2,3	Hp	YSI 400 Serie	6	
	Siemens, Dräger, Hellige	PH-AR-275A	Adult	Rectal/Esophageal	x		275	4,0	Phone	YSI 400 Serie	4,8	
		PH-PR-275A	Pediatric	Rectal/Esophageal	x		275	3,2	Phone	YSI 400 Serie	4	
		PH-NR-275A	Infant/Neonatal	Rectal/Esophageal	x		275	2,3	Phone	YSI 400 Serie	2,7	
	Siemens, Dräger, Hellige etc.	PH-SAR-275A	Adult	Skin-Contact	x		275	4,0	Phone	YSI 400 Serie	18	
		PH-SPR-275A	Pediatric	Skin-Contact	x		275	3,2	Phone	YSI 400 Serie	10	
		PH-SNR-275A	Infant/Neonatal	Skin-Contact	x		275	2,3	Phone	YSI 400 Serie	6	
	Siemens, Dräger, Hellige	PHR-AR-275A	Adult	Rectal/Esophageal	x		275	4,0	Phone right-angled	YSI 400 Serie	5	
		PHR-PR-275A	Pediatric	Rectal/Esophageal	x		275	3,2	Phone right-angled	YSI 400 Serie	4	
		PHR-NR-275A	Infant/Neonatal	Rectal/Esophageal	x		275	2,3	Phone right-angled	YSI 400 Serie	3	
	Siemens, Dräger, Hellige etc.	PHR-SAR-275A	Adult	Skin-Contact	x		275	4,0	Phone right-angled	YSI 400 Serie	18	
		PHR-SPR-275A	Pediatric	Skin-Contact	x		275	3,2	Phone right-angled	YSI 400 Serie	10	
		PHR-SNR-275A	Infant/Neonatal	Skin-Contact	x		275	2,3	Phone right-angled	YSI 400 Serie	6	
	Siemens SC 6000, SC 9000	SI-AR-275A	Adult	Rectal/Esophageal	x		275	4,0	Siemens-Multimed	YSI 400 Serie	4,8	
		SI-PR-275A	Pediatric	Rectal/Esophageal	x		275	3,2	Siemens-Multimed	YSI 400 Serie	4	
		SI-NR-275A	Infant/Neonatal	Rectal/Esophageal	x		275	2,3	Siemens-Multimed	YSI 400 Serie	2,7	
	Siemens SC 6000, SC 9000	SI-SAR-275A	Adult	Skin-Contact	x		275	4,0	Siemens-Multimed	YSI 400 Serie	18	
		SI-SPR-275A	Pediatric	Skin-Contact	x		275	3,2	Siemens-Multimed	YSI 400 Serie	10	
		SI-SNR-275A	Infant/Neonatal	Skin-Contact	x		275	2,3	Siemens-Multimed	YSI 400 Serie	6	